CleanAIR®

CleanAIR®
PERSONAL RESPIRATORY PROTECTION SYSTEMS

USER MANUAL

CleanAIR® MedicAER®

Malina Safety®

MALINA - Safety s.r.o.
Luční 11,
466 01 Jablonec n. Nisou
Czech Republic


Tel. +420 483 356 600
export@malina-safety.cz
www.malina-safety.com

ENG	4	DUT	31	GER	59	NOR	86	SWE	114
CZE	11	SPA	38	HUN	66	POL	93	TUK	121
CHI	18	FIN	45	HEB	73	POR	100		
DAN	24	FRE	52	ITA	79	RUS	107		

NA-068-R02

Possible headtop combinations

CA-1	CA-2	CA-10
TH3	TH3	TH3

UniMask	CA-3 acetate visor	
TH3	TH2	

Shigematsu GX02	Shigematsu CF02	CA-5
TM3	TM3	TM2

**CleanAIR®****MedicAER®**

Powered Air Purifying Respirator (PAPR)

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Contents:

1. Introduction
2. Instructions for use
3. Unpacking and assembly
4. Before use
5. Maintenance, cleaning and decontamination
6. Spare parts and their replacement
7. Possible faults
8. Basic description of unit control and settings
9. Storage
10. Warranty
11. Technical data
12. List of spare parts / Accessories

1. Introduction

CleanAIR® is a personal respiratory protection system based on overpressure of filtered air in the breathing zone. The powered air purifying respirator (hereinafter referred to as the "PAPR") sucks in air from the surrounding environment and supplies it to the protection hood or mask. The resulting overpressure prevents contaminants from passing into the breathing zone. At the same time the light overpressure ensures high user comfort even when worn for a long time, because there is no need to overcome any breathing resistance of the filter. The system keeps the selected airflow regardless of the filter clogging or the degree of battery charge.

To ensure efficient protection, it is necessary to choose a suitable combination of PAPR unit and protective headtop. It is also important to use the proper filters of the same type in regards to the kind of contamination and its concentration.

The CleanAIR® MedicAER® meets the requirements of the EN 12941 and EN 12942. Depending on the filters used, it provides protection against non-toxic as well as toxic particles and harmful gases.

2. Instructions for use

Carefully read and follow the instructions stated in this user manual. The user must know exactly how to use the protective equipment correctly.

- If, for any reason, the PAPR unit stops working, the user must leave the contaminated environment immediately. It applies even though the system provides sufficient protection according to EN 12942 while switched off.
- If used with the hood when the unit is switched off, the respiratory protection is low or none. In that case there is a risk of rapid increase of concentration of carbon dioxide and decrease of oxygen concentration. The situation when the unit is switched off when used with the hood is abnormal.
- CleanAIR® MedicAER® must not be used when the amount of incoming air is insufficient.
- CleanAIR® MedicAER® must not be used in areas where oxygen concentration is below 17%.
- CleanAIR® MedicAER® must not be used in areas where the user has no knowledge about the nature of contamination or its concentration.
- CleanAIR® MedicAER® must not be used in areas where there is a risk of explosion.
- CleanAIR® MedicAER® must not be used in confined spaces, such as closed tanks, tunnels, sewers, etc.
- Before each use, check the airflow with the enclosed flow indicator. The flow rate is insufficient if the top edge of the cone is in the red zone (see Section 4.2 "Airflow test").
- Worn headtops must fit tightly onto face to ensure declared level of protection for the wearer. If a beard or long hair interferes with the sealing line, this may negatively affect the protection provided by the system.

 It is very important to choose appropriate filters for the relevant kind of contamination.
Follow the instructions in the user manual enclosed with the filters.

- The unit provides protection against different contaminants depending on selected filters.
- Filters designed for filtering solid and liquid aerosols do not protect the user against gases.
- Filters designed for filtering gases and vapours do not protect the user against particles.
- Filters designed just for filtering gases and vapours are not recommended to use on PAPR unit if they are not combined with a particle insert for the filtration of particles.
- Filters must not be attached directly to the hood.
- Use only original CleanAIR® particle or combined filters designed for your CleanAIR® MedicAER® PAPR unit.
- Always replace the filters if you notice a change in the smell of the air coming out of the unit.
- CleanAIR® MedicAER® does not remove carbon monoxide (CO) or carbon dioxide (CO₂) from the air!
- Filters designed for filtering mercury are prohibited to use for longer than 50 working hours.
- Always follow your local regulations for using PAPR units, filters and PPE.

 If any rule or instruction stated in this manual is not observed, the warranty becomes null and void!

3. Unpacking / Assembly

3.1. Unpacking the CleanAIR® MedicAER®

Check whether the package contents are complete and that nothing has been damaged during transport. The complete system, including accessories, contains the following components:

	51 MO 00FC
MedicAER® unit	✓
Exchangeable battery	✓
Comfort padded belt	✓
Light flexi hose	✓
Case	✓
Flow indicator	✓
Battery charger	✓
User manual	✓

3.2. Assembly

1. Insert the battery into the unit.
2. Connect the hose to the PAPR unit and make sure the connection is tight enough.
3. Attach filters to the unit - make sure they are both the same type!
4. Connect the hose to the headtop.

4. Before use

4.1. Control procedure before every use - make sure that:

- All components are in good condition without any visible damage. Replace every damaged and worn out part.
- The hose is connected correctly to the PAPR unit as well as to the headtop.
- Airflow is sufficient (Section 4.2.).

4.2. Airflow test

1. Disconnect the air hose from the unit.
2. Connect the flow indicator to the unit.
3. Switch the PAPR unit on.

 If the pointer's top edge enters the red zone, airflow is insufficient and the filters must be replaced!

4.3. Alarm test

To check if the alarm works properly, close the air output with your palm and turn the unit on. You should hear and see alarm signals within 10 seconds.

5. Maintenance, cleaning and decontamination

At the end of each working shift, it is recommended to check and clean all parts of the system and replace those that are damaged.

- Do not use aggressive detergents and solvents for cleaning! We recommend using common non-abrasive cleaning products.
- Cleaning products must never penetrate the PAPR unit or the battery!
- Use a damp cloth for cleaning and wipe the surface dry afterwards.
- Cleaning should be done in a well ventilated room.

 Beware of breathing in any particles and fibres settled on parts of the PAPR unit or its accessories!

6. Spare parts and their replacement

6.1. Replacing filters

Only filters designed for CleanAIR® MedicAER® may be used! The filters must be new, unused and in the original packaging. The expiration date must not have elapsed – it is marked on a sticker on the body of the filter.

Filter replacement:

1. Unscrew the filter from the unit.
2. Check the rubber sealing ring at the inputs to the unit for possible damage.
3. Fasten the new filters onto the unit and screw them tight.

There are two filter threads RD40x1/7" on the unit [standardized according to EN 148-1]. Whenever changing filters, you should always change a set of two filters of the same type at the same time!

It is possible to use a pre-filter. This catches larger particles and especially aerosols, which can clog up the filter during paint-spraying. Thus a pre-filter prolongs the life-time of the filter.

Install only new filters without any visible damage.

It is forbidden to clean a filter or blow through it in any way!

From the hygiene point of view it is not recommended to use filters for more than 1 month after their first use.

6.2. Batteries

For technical reasons, batteries are not fully charged when leaving the manufacturer. They are supplied with a 50 - 70% charge. Therefore either fully charge the battery before first use, or expect shorter runtime in the first cycle. For maximum battery life, follow the storage instructions (see Section 9. "Storage")

6.2.1. Battery charging

- Slide the battery out of the unit.
- Connect the charger to the 100 - 230 V ~ 50/60 Hz mains.
- Connect the battery to the charger. Charging is indicated by a red light, the green light indicates connection to power. Charging time is < 3 hours.
- When the red light goes off, the battery is fully charged.
- After charging, disconnect the battery from the charger and the charger from the mains.

It is also possible to charge the battery while the battery is still plugged into the unit. The unit must be turned off during charging.

It is not recommended to leave the charger connected to the mains when not in use.

Charging starts after connecting the battery using the charger to the power supply. After the battery is fully charged, the charger switches to maintenance mode to keep the battery fully charged. The charging time is less than 3 hours.

Red LED on: fast charging
Red LED off: charge complete

Red LED flashing while plugged: error
Red LED flashing while unplugged: pre-charge detection

The battery charger is for indoor use only. Never charge the battery in potentially explosive atmosphere.
It is forbidden to use the battery charger for other purposes than those determined by the producer.

6.2.2. Battery replacement

see pictorial annex "How to insert the battery" and "How to remove the battery"

1. Grip the powered respiratory unit by the rim on the top side of the battery with your left hand. Release the safety lock located on the bottom of the unit between the filters. Then push the battery out of the unit using your thumb (press the rubber place on the battery).
2. Insert the battery into the corresponding opening of the unit until the click-in stop locks the battery in the right position.


7. Possible faults and alarms

In case of any fault, such as a sudden decrease or increase in air supply when the user is in a contaminated working area, it is necessary to leave the worksite immediately.


If there is a problem, check the following:

- Whether the unit is assembled correctly.
- Condition of the battery.
- Filter clogging.
- Whether the sealing of the safety hood provides a good fit.


The unit is equipped with warning systems that are activated in the following situations:

Blocked filter
Change filter

Low battery
Charge battery

Time expiration
Change filter

Blocked filter - if the unit cannot maintain the chosen airflow, the electronics will automatically decrease the airflow by one level and you will hear an acoustic signal. If the unit is unable to maintain the lowest possible airflow, you will be warned by an audiovisual alarm.

Low battery - an alarm notifies the user when the battery has run out. It is best to leave the contaminated area and replace or recharge the battery.




Filter expiration – the alarm should be pre-set manually before you start using new filters. Once you start using the filters, the unit starts the countdown from the pre-set time. When the pre-set time expires, the unit will start to warn you.

8. Basic description of unit control and settings

Switching on

After turning the unit on, you will see information about the manufacturer of the unit on your display.

Main screen information

The next screen shows current airflow, filter clogging and battery level. You can also see the current headtop and filter setting.

Airflow adjustment

You can adjust airflow by simple pressing the \oplus or \ominus button on the side of the unit. The available airflow setting depends on the headtop and filter setting.

Rotating display

By briefly pressing down on both buttons you can rotate the display clockwise by 90 degrees.

Standard
Garment
Filter
Particle
Filter Timer
Alarm Off
Reser timer
No
> Language <
English
Display View
Full only
Default settings
No

How to get to the main menu

Press the \oplus and \ominus button at the same time and hold them for 2 seconds. The main menu will appear.

How to navigate in the menu and change settings:

- Press the \ominus button and jump to the next menu item, or hold the \ominus button for 2 sec to jump to the previous item.
- Press the \oplus button and change the settings of the item.
- Press both buttons for 2 seconds to save the settings and leave the menu [not pressing anything for 5 sec does the same].

Unit headtop setting:

Go to the unit menu (\oplus and \ominus for 2 sec) and repeat pressing \oplus until the required headtop standard appears.

9. Storage

All parts of the CleanAIR® system must be stored at temperatures between -10 °C and 55 °C, with relative air humidity between 20 and 95% Rh.

Batteries will self-discharge during storage. Therefore, it is recommended to charge the battery for 1 hour at least every 3 months. The optimal battery charge level for longer storage is between 50 - 70% of capacity. After longer storage time it is necessary to repeat 3 charging cycles to reach battery's full capacity.

10. Warranty

The warranty ensures that you will receive a replacement if a product has any manufacturing or material defects that appear within 12 months of the date of purchase. The same period for batteries is 6 months from the date of purchase. The warranty claim must be reported to the sales department/retailer. At the same time, proof of purchase must be submitted (i.e. an invoice or certificate of delivery). The warranty can only be acknowledged if no interventions into the powered respiratory unit and charger have been made. The warranty does not cover defects caused by the late replacement of filters, or by using a filter which was damaged by cleaning or being blown through.

The manufacturer strongly recommends to perform periodical inspections in CleanAIR® authorized service center. Inspections should be done according the local law and at least once in every two years.

11. Technical data

Technical specifications - CleanAIR® MedicAER®

Airflow	120 - 235 l/min
Display	TFT display
Weight of unit, including battery	960 g
Unit noise	55 to 61 dB
Battery lifespan	Max. 500 recharging cycles
One charging cycle	less than 3 hours
Belt - waist size	up to 1,500 mm
Operating temperature	0 to 60 °C
Operating humidity	20 to 95% Rh
Certification	EN 12941/A2 TH2/TH3; EN 12942/A2 TM2/TM3

The unit conforms with IP65 - IP68 (Ingress Protection against liquids and solid particles):

- IP64: When the unit is powered on [no additional requirements]
- IP65: When the unit is on [attached spark arrestor and pre-filter holder set - 50 0112]
- IP68: When unit is off [with inlets and outlets closed by decontamination plugs - 51 00 46]

Notified body for testing CE:

Occupations Safety Research Institute
Jeruzálemská 1283/9, 110 00 Praha 1
Autorizovaná osoba 235
Notified body 1024

Declaration of Conformity is available at: <https://www.clean-air.cz/doc/MedicAER.pdf>

12. List of spare parts

Order code:	Description:
51 MO 00FC	CleanAIR® MedicAER® [with charger, comfort belt, standard battery, air hose, set of P R SL filters, flow indicator and plastic transport case]
51 MO 10	Standard battery for MedicAER® [Li-Ion 14,4 V / 2,6 Ah]
51 MO 20	Heavy Duty battery for MedicAER® [Li-Ion 14,4 V / 5,2 Ah]
52 00 44.1	Comfort harness
51 00 41	Decontaminable belt
51 00 42	Decontaminable harness
51 00 30EUR	Charger Li-Ion 14,4 V - EURO plug
51 00 30UK	Charger Li-Ion 14,4 V - UK plug
71 00 60	Light flexi hose QuickLOCK
71 00 60L	Light flexi hose QuickLOCK extended
71 00 86	Rubber hose QuickLOCK
71 00 86L	Rubber hose QuickLOCK extended
70 00 62F	Hose cover flameproof
70 00 62C	Hose cover chemical resistant
71 00 92	Comfort padded belt
52 00 90	Flow indicator CA40x1/7"
51 00 46	Set of decontamination plugs

CleanAIR® MedicAER® - compatible filters

Order number	Filter type *	Colour code	Main application
50 00 48	P3	White	Particles in the form of solid and liquid aerosols (dust, fibres, toxic and non-toxic fumes, bacteria and viruses)
50 03 57	A1P3	Brown White	Organic gases and vapours with a boiling point above 65 °C Particles in the form of solid and liquid aerosols
50 03 64	A1B1E1P3	Brown Grey Yellow White	Organic gases and vapours with a boiling point above 65 °C Inorganic gases and vapours Acid gases and vapours Particles in the form of solid and liquid aerosols
50 01 57	A2P3	Brown White	Organic gases and vapours with a boiling point above 65 °C Particles in the form of solid and liquid aerosols
50 01 67	A2B2P3	Brown Grey White	Organic gases and vapours with a boiling point above 65 °C Inorganic gases and vapours Particles in the form of solid and liquid aerosols
50 01 68	A2B2E2K2P3	Brown Grey Yellow Green White	Organic gases and vapours with a boiling point above 65 °C Inorganic gases and vapours Acid gases and vapours Ammonia and organic amines Particles in the form of solid and liquid aerosols
50 01 66	A2B2E2K2HgP3	Brown Grey Yellow Green Red White	Organic gases and vapours with a boiling point above 65 °C Inorganic gases and vapours Acid gases and vapours Ammonia and organic amines Mercury Particles in the form of solid and liquid aerosols

* Filters for CleanAIR® are labelled according to standards EN 14387:2004 + A1, EN 143:2000 / A1:2006, EN 12941/A2 and EN 12942/A2. This table does not contain a complete list of compatible filters. For a complete list of filters, please contact your supplier for more information.

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